

For Research Use Only

PPP2R3C Recombinant antibody, PBS Only

Catalog Number: 86449-1-PBS



Basic Information

Catalog Number:

86449-1-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19282

GenBank Accession Number:

BC063438

GeneID (NCBI):

55012

UNIPROT ID:

Q969Q6

Full Name:

protein phosphatase 2 (formerly 2A),
regulatory subunit B", gamma

Calculated MW:

453 aa, 53 kDa

Observed MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

251047G3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86449-1-RR (PPP2R3C antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86449-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86449-1-RR against Human PPP2R3C were performed. The affinity constant is below 1 pM.