

For Research Use Only

BRSK1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86462-2-PBS



Basic Information

Catalog Number:

86462-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3357

GenBank Accession Number:

BC016681

GeneID (NCBI):

84446

UNIPROT ID:

Q8TDC3

Full Name:

BR serine/threonine kinase 1

Calculated MW:

87 kDa

Purification Method:

Protein A purification

CloneNo.:

251341E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BRSK1 concentration was determined to be 8.4 ng/mL in Mouse brain tissue extract based on a 7.0 mg/mL extract load.



Cytometric bead array standard curve of MP02463-1, BRSK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86462-2-PBS. Detection antibody: 86462-1-PBS. Standard: Ag3357. Range: 0.312-40 ng/mL.