

For Research Use Only

AGPS Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86472-2-PBS



Basic Information

Catalog Number: 86472-2-PBS	GenBank Accession Number: BC141820	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 8540	CloneNo.: 251340E1
Isotype: IgG	UNIPROT ID: O00116	
Immunogen Catalog Number: AG15233	Full Name: alkylglycerone phosphate synthase	
	Calculated MW: 658 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean AGPS concentration was determined to be 110.4 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02459-1, AGPS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86472-2-PBS. Detection antibody: 86472-1-PBS. Standard: Ag15233. Range: 0.625-40 ng/mL.