

For Research Use Only

BAF170 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86473-1-PBS



Basic Information

Catalog Number:

86473-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2634

GenBank Accession Number:

BC013045

GeneID (NCBI):

6601

UNIPROT ID:

Q8TAQ2

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily c, member 2

Calculated MW:

1214 aa, 132 kDa

Purification Method:

Protein A purification

CloneNo.:

251330B10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BAF 170 concentration was determined to be 7.4 ng/mL in HepG2 cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP02458-1, BAF170 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86473-2-PBS. Detection antibody: 86473-1-PBS. Standard: Ag2634. Range: 0.625-80 ng/mL.