

For Research Use Only

# BAF170 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86473-4-PBS



## Basic Information

Catalog Number:

86473-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2634

GenBank Accession Number:

BC013045

GeneID (NCBI):

6601

UNIPROT ID:

Q8TAQ2

Full Name:

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily c, member 2

Calculated MW:

1214 aa, 132 kDa

Purification Method:

Protein A purification

CloneNo.:

251330C4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean BAF 170 concentration was determined to be 657.43 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load and 668.02 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02458-2, Human BAF 170 Recombinant Matched Antibody Pair - PBS only. 86473-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2634. 86473-2-PBS was HRP conjugated as the detection antibody. Range: 6.25-200 ng/mL