

For Research Use Only

CHMP6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86474-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
86474-1-PBS	BC010108	Protein A purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	79643	251259B2
Source:	UNIPROT ID:	
Rabbit	Q96FZ7	
Isotype:	Full Name:	
IgG	chromatin modifying protein 6	
Immunogen Catalog Number:	Calculated MW:	
AG9357	201 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CHMP6 concentration was determined to be 220.3 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02462-1, CHMP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86474-2-PBS. Detection antibody: 86474-1-PBS. Standard: Ag9357. Range: 0.312-40 ng/mL.