

For Research Use Only

CHMP6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86474-2-PBS



Basic Information

Catalog Number:

86474-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9357

GenBank Accession Number:

BC010108

GeneID (NCBI):

79643

UNIPROT ID:

Q96FZ7

Full Name:

chromatin modifying protein 6

Calculated MW:

201 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

251259D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CHMP6 concentration was determined to be 220.3 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02462-1, CHMP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86474-2-PBS. Detection antibody: 86474-1-PBS. Standard: Ag9357. Range: 0.312-40 ng/mL.