

For Research Use Only

UMP/CMP kinase Recombinant monoclonal antibody, PBS Only

Catalog Number: 86475-3-PBS



Basic Information

Catalog Number: 86475-3-PBS	GenBank Accession Number: BC014961	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 51727	CloneNo.: 251233A1
Isotype: IgG	UNIPROT ID: P30085	
Immunogen Catalog Number: AG1920	Full Name: cytidine monophosphate (UMP-CMP) kinase 1, cytosolic	
	Calculated MW: 228 aa, 26 kDa	
	Observed MW: 22 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

Pyrimidine nucleoside monophosphate kinase [UMP/CMP kinase (UMP/CMK)] (CMK1) is also named as CMK, CMK1, UCK, UCK1, UCK2, UCK3 and belongs to the adenylate kinase family. It plays a crucial role in the formation of UDP, CDP, and dCDP, which are required for cellular nucleic acid synthesis (PMID:11912132). CMK1 can use a broad spectrum of phosphate donors but the best phosphate donors are ATP and dATP. CMK1 has been reported as a prognostic marker for non-small cell lung cancer, pancreatic cancer and breast cancer (PMID: 21642870, 22838950, 27558661). The longest deduced protein contains 228 amino acids and has a calculated molecular mass of 26 kD. However, translation likely begins at the second initiation methionine, resulting in a protein of 196 amino acids (PMID:10462544). The shorter form is the endogenously translated protein (PMID:11912132). CMK1 has 3 isoforms with the molecular mass of 22, 16 and 26 kDa.

Storage

Storage:
Store at -80°C.


The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86475-3-RR (CMPK1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86475-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86475-3-RR against Human UMP/CMP kinase were performed. The affinity constant is below 1 pM.