

For Research Use Only

CD20 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86480-2-PBS



Basic Information

Catalog Number:

86480-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17309

GenBank Accession Number:

BC002807

GeneID (NCBI):

931

UNIPROT ID:

P11836-1

Full Name:

membrane-spanning 4-domains,
subfamily A, member 1

Calculated MW:

297 aa, 33 kDa

Purification Method:

Protein A purification

CloneNo.:

260788G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD20 concentration was determined to be 20.5 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load, HepG2 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03555-1, CD20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86480-10-PBS. Detection antibody: 86480-2-PBS. Standard: Ag17309. Range: 0.469-60 ng/mL.