

For Research Use Only

CREBBP Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86481-2-PBS



Basic Information

Catalog Number: 86481-2-PBS	GenBank Accession Number: NM_001079846	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 1387	CloneNo.: 253020G10
Isotype: IgG	UNIPROT ID: Q92793	
Immunogen Catalog Number: AG27489	Full Name: CREB binding protein	
	Calculated MW: 265 aa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CREBBP concentration was determined to be 17.6 ng/mL in HEK-293T cell extract based on a 1.2 mg/mL extract load, 11.2 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load and 2.8 ng/mL in PC-3 cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP03100-2, Human CREBBP Recombinant Matched Antibody Pair - PBS only. 86481-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27489. 86481-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.