

For Research Use Only

CD146/MCAM Recombinant monoclonal antibody, PBS Only

Catalog Number: 86492-1-PBS



Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 86492-1-PBS | GenBank Accession Number: NM_006500.2 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 4162 | CloneNo.: 250530B11 |
| Isotype: IgG | UNIPROT ID: P43121-1 | |
| Immunogen Catalog Number: EG2938 | Full Name: melanoma cell adhesion molecule | |
| | Calculated MW: 72 kDa | |
| | Observed MW: 120 kDa | |

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

CD146, also known as melanoma cell adhesion molecule (MCAM) or MUC18, originally identified as a biomarker of melanoma progression, is a transmembrane glycoprotein of 113-130 kDa, belonging to the immunoglobulin (Ig) superfamily (PMID: 8378324; 25993332). Structurally, it consists of five Ig domains, a transmembrane domain, and a cytoplasmic region. In normal adult tissue, CD146 is primarily expressed by vascular endothelium and smooth muscle. CD146 is a key cell adhesion protein in vascular endothelial cell activity and angiogenesis, and has been used as marker of circulating endothelial cells (CECs) (PMID: 19356677). In addition to the membrane-anchored form of CD146, a soluble form of CD146 (sCD146, 105 kDa) has also been found in human plasma and in the supernatant of cultured human endothelial cells (PMID: 9462829; 19229070; 16374253; 14597988).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86492-1-RR (MCAM antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86492-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HUVEC cells using CD146/MCAM antibody (86492-1-RR, Clone: 250530B11) at dilution of 1:800 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86492-1-PBS in a different storage buffer formulation.