

For Research Use Only

GTF2F1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86494-1-PBS



Basic Information

Catalog Number:

86494-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0133

GenBank Accession Number:

BC000120

GeneID (NCBI):

2962

UNIPROT ID:

P35269

Full Name:

general transcription factor IIF,
polypeptide 1, 74kDa

Calculated MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

251305C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GTF2F1 concentration was determined to be 28.6 ng/mL in Mouse lung tissue extract based on a 5.9 mg/mL extract load.



Cytometric bead array standard curve of MP02475-1, GTF2F1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86494-2-PBS. Detection antibody: 86494-1-PBS. Standard: Ag0133. Range: 0.625-80 ng/mL.