

For Research Use Only

# TPOR/MPL Recombinant monoclonal antibody, PBS Only

Catalog Number:86495-3-PBS



## Basic Information

Catalog Number:

86495-3-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC153092

GeneID (NCBI):

4352

UNIPROT ID:

P40238

Full Name:

myeloproliferative leukemia virus oncogene

Calculated MW:

635 aa, 71 kDa

Observed MW:

84 kDa, 74 kDa

Purification Method:

Protein A purification

CloneNo.:

251731G4

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



HEL cells were subjected to SDS PAGE followed by western blot with 86495-3-RR (MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86495-3-PBS in a different storage buffer formulation.



human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 86495-3-RR (MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86495-3-PBS in a different storage buffer formulation.