

For Research Use Only

TPOR/MPL Recombinant monoclonal antibody

Catalog Number:86495-3-RR



Basic Information

Catalog Number:

86495-3-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC153092

GeneID (NCBI):

4352

UNIPROT ID:

P40238

Full Name:

myeloproliferative leukemia virus oncogene

Calculated MW:

635 aa, 71 kDa

Observed MW:

84 kDa, 74 kDa

Purification Method:

Protein A purification

CloneNo.:

251731G4

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human peripheral blood platelets, HEL cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEL cells were subjected to SDS PAGE followed by western blot with 86495-3-RR (MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 86495-3-RR (MPL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.