

For Research Use Only

# USP16 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86498-2-PBS



## Basic Information

<b>Catalog Number:</b> 86498-2-PBS	<b>GenBank Accession Number:</b> BC030777	<b>Purification Method:</b> Protein A purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 10600	<b>CloneNo.:</b> 251272E8
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9Y5T5	
<b>Isotype:</b> IgG	<b>Full Name:</b> ubiquitin specific peptidase 16	
<b>Immunogen Catalog Number:</b> AG4981	<b>Calculated MW:</b> 94 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP02643-1, Human USP16 Recombinant Matched Antibody Pair - PBS only. 86498-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag4981. 86498-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean USP16 concentration was determined to be 24.30 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load and 64.16 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.