

For Research Use Only

USP16 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86498-3-PBS



Basic Information

Catalog Number: 86498-3-PBS	GenBank Accession Number: BC030777	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 10600	CloneNo.: 251272E4
Source: Rabbit	UNIPROT ID: Q9Y5T5	
Isotype: IgG	Full Name: ubiquitin specific peptidase 16	
Immunogen Catalog Number: AG4981	Calculated MW: 94 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02643-1, Human USP16 Recombinant Matched Antibody Pair - PBS only. 86498-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag4981. 86498-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean USP16 concentration was determined to be 24.30 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load and 64.16 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.