

For Research Use Only

FABP3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86511-1-PBS



Basic Information

Catalog Number: 86511-1-PBS	GenBank Accession Number: BC007021	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 2170	CloneNo.: 251339B6
Source: Rabbit	UNIPROT ID: P05413	
Isotype: IgG	Full Name: fatty acid binding protein 3, muscle and heart (mammary-derived growth inhibitor)	
Immunogen Catalog Number: AG1069	Calculated MW: 15 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FABP3 concentration was determined to be 29.2 ug/mL in Human milk.



Cytometric bead array standard curve of MP02478-1, FABP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86511-2-PBS. Detection antibody: 86511-1-PBS. Standard: Ag1069. Range: 1.25-80 ng/mL