

For Research Use Only

KHDRBS3 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86516-2-PBS



Basic Information

Catalog Number: 86516-2-PBS	GenBank Accession Number: BC032606	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10656	CloneNo.: 251366G10
Isotype: IgG	UNIPROT ID: O75525	
Immunogen Catalog Number: AG4438	Full Name: KH domain containing, RNA binding, signal transduction associated 3	
	Calculated MW: 346 aa, 38, 8 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean KHDRBS3 concentration was determined to be 756.4 ng/mL in Mouse brain tissue extract based on a 7.0 mg/mL extract load.



Cytometric bead array standard curve of MP02486-1, KHDRBS3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86516-2-PBS. Detection antibody: 86516-1-PBS. Standard: Ag4438. Range: 0.781-100 ng/mL.