

For Research Use Only

ALDH1A2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86518-3-PBS



Basic Information

Catalog Number:

86518-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5062

GenBank Accession Number:

BC030589

GeneID (NCBI):

8854

UNIPROT ID:

O94788

Full Name:

aldehyde dehydrogenase 1 family, member A2

Calculated MW:

518 aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

251371E6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ALDH1A2 concentration was determined to be 2,903.03 ng/mL in SKOV-3 cell extract based on a 1.1 mg/mL extract load.



Sandwich ELISA standard curve of MP02658-1, Human ALDH1A2 Recombinant Matched Antibody Pair - PBS only. 86518-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5062. 86518-3-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL