

For Research Use Only

SH3D19 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86522-1-PBS



Basic Information

Catalog Number: 86522-1-PBS	GenBank Accession Number: BC108890	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 152503	CloneNo.: 251364F3
Isotype: IgG	UNIPROT ID: Q5HYK7	
Immunogen Catalog Number: AG16200	Full Name: SH3 domain containing 19	
	Calculated MW: 790 aa, 87 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean SH3D19 concentration was determined to be 1.4 ng/mL in Y79 cell extract based on a 1.1 mg/mL extract load, 1.3 ng/mL in HT-1080 cell extract based on a 1.4 mg/mL extract load.



Cytometric bead array standard curve of MP02487-1, SH3D19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86522-2-PBS. Detection antibody: 86522-1-PBS. Standard: Ag16200. Range: 0.312-40 ng/mL.