

For Research Use Only

# REXO4 Recombinant antibody, PBS Only

Catalog Number: 86528-1-PBS



## Basic Information

Catalog Number:

86528-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13473

GenBank Accession Number:

BC009274

GeneID (NCBI):

57109

UNIPROT ID:

Q9GZR2

Full Name:

REX4, RNA exonuclease 4 homolog (S. cerevisiae)

Calculated MW:

422 aa, 47 kDa

Observed MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

251440H2

## Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86528-1-RR (REXO4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86528-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using REXO4 antibody (86528-1-RR, Clone: 251440H2) at dilution of 1:1500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 86528-1-PBS in a different storage buffer formulation.