

For Research Use Only

Shootin 1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86541-1-PBS



Basic Information

Catalog Number: 86541-1-PBS	GenBank Accession Number: BC022348	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 57698	CloneNo.: 251435A2
Source: Rabbit	UNIPROT ID: A0MZ66	
Isotype: IgG	Full Name: KIAA1598	
Immunogen Catalog Number: AG18294	Calculated MW: 72 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Shootin 1 concentration was determined to be 50.1 ng/mL in SH-SY5Y cell extract based on a 1.4 mg/mL extract load, 47.7 ng/mL in Neuro-2a cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02499-1, Shootin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86541-2-PBS. Detection antibody: 86541-1-PBS. Standard: Ag18294. Range: 0.781-100 ng/mL.