

For Research Use Only

KYNU Recombinant antibody

Catalog Number: 86555-3-RR



Basic Information

Catalog Number:

86555-3-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2365

GenBank Accession Number:

BC000879

GeneID (NCBI):

8942

UNIPROT ID:

Q16719

Full Name:

kynureninase (L-kynurenine hydrolase)

Calculated MW:

465aa, 52 kDa; 307aa, 35 kDa

Observed MW:

45-52 kDa

Purification Method:

Protein A purification

CloneNo.:

251437D1

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, A549 cells

Background Information

KYNU(kynureninase) is also named as L-kynurenine hydrolase and belongs to the kynureninase family. It is a pyridoxal-5'-phosphate-(pyridoxal-P)-dependent enzyme that catalyses the cleavage of L-kynurenine and L-3-hydroxykynurenine into anthranilic and 3-hydroxyanthranilic acids, respectively(PMID:8706755). In SDS-PAGE, it can detect a 48 kDa band lower than the predicted molecular weight of 52 kDa due to limited proteolysis of the enzyme which might occur during handling of the protein(PMID:8706755).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86555-3-RR (KYNu antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.