

For Research Use Only

RCAN2 Recombinant antibody, PBS Only

Catalog Number: 86561-1-PBS



Basic Information

Catalog Number:

86561-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3597

GenBank Accession Number:

BC038509

GeneID (NCBI):

10231

UNIPROT ID:

Q14206

Full Name:

regulator of calcineurin 2

Calculated MW:

225 aa, 25 kDa

Observed MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

251381D10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86561-1-RR (RCAN2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86561-1-PBS in a different storage buffer formulation.