

For Research Use Only

RCAN2 Recombinant antibody

Catalog Number: 86561-1-RR



Basic Information

Catalog Number:

86561-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3597

GenBank Accession Number:

BC038509

GeneID (NCBI):

10231

UNIPROT ID:

Q14206

Full Name:

regulator of calcineurin 2

Calculated MW:

225 aa, 25 kDa

Observed MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

251381D10

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat, pig

Positive Controls:

WB : mouse brain tissue, rat brain tissue, pig brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86561-1-RR (RCAN2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.