

For Research Use Only

RCAN2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86561-2-PBS



Basic Information

Catalog Number: 86561-2-PBS	GenBank Accession Number: BC038509	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10231	CloneNo.: 251381F2
Isotype: IgG	UNIPROT ID: Q14206	
Immunogen Catalog Number: AG3597	Full Name: regulator of calcineurin 2	
	Calculated MW: 225 aa, 25 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of fourteen individual healthy human donors was measured. The human RCAN2 concentration of detected samples was determined to be 6.78 ng/mL with a range of 0.58 - 29.64 ng/mL.



Sandwich ELISA standard curve of MP02659-1, Human RCAN2 Recombinant Matched Antibody Pair - PBS only. 86561-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3597. 86561-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL