

For Research Use Only

Mouse IL-17C Recombinant antibody, PBS Only

Catalog Number: 86563-1-PBS



Basic Information

Catalog Number:

86563-1-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2840

GenBank Accession Number:

NM_145834.3

GeneID (NCBI):

234836

UNIPROT ID:

Q8K4C5

Full Name:

interleukin 17C

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

250601E2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with 86563-1-RR (IL17c antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86563-1-PBS in a different storage buffer formulation.