

For Research Use Only

Mouse Il17c Recombinant antibody, PBS Only (Capture)

Catalog Number:86563-3-PBS



Basic Information

Catalog Number: 86563-3-PBS	GenBank Accession Number: NM_145834.3	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 234836	CloneNo.: 250601G8
Source: Rabbit	UNIPROT ID: Q8K4C5	
Isotype: IgG	Full Name: interleukin 17C	
Immunogen Catalog Number: EG2840	Calculated MW: 22 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IL-17C concentration was determined to be 3.61 ng/mL in mouse intestine tissue extract based on a 4.8 mg/mL extract load.



Sandwich ELISA standard curve of MP02646-1, Mouse IL-17C Recombinant Matched Antibody Pair - PBS only. 86563-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2840. 86563-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL