

For Research Use Only

PCDH17 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86566-3-PBS



Basic Information

Catalog Number: 86566-3-PBS	GenBank Accession Number: BC028165	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 27253	CloneNo.: 251404C10
Source: Rabbit	UNIPROT ID: O14917	
Isotype: IgG	Full Name: protocadherin 17	
Immunogen Catalog Number: EG4889	Calculated MW: 1159 aa, 126 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PCDH17 concentration was determined to be 20.8 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02677-1, Human PCDH17 Recombinant Matched Antibody Pair - PBS only. 86566-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4889. 86566-3-PBS was HRP conjugated as the detection antibody. Range: 97.7-6250 pg/mL.