

For Research Use Only

# ZDHHC8 Recombinant monoclonal antibody, PBS Only

Catalog Number:86572-2-PBS



## Basic Information

Catalog Number:

86572-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26028

GenBank Accession Number:

BC053544

GeneID (NCBI):

29801

UNIPROT ID:

Q9ULC8

Full Name:

zinc finger, DHHC-type containing 8

Calculated MW:

83 kDa

Observed MW:

70-85 kDa

Purification Method:

Protein A purification

CloneNo.:

251438A7

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86572-2-RR (ZDHHC8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86572-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86572-2-RR against Human ZDHHC8 were performed. The affinity constant is 0.915 nM.