

For Research Use Only

ZDHHC8 Recombinant monoclonal antibody

Catalog Number: 86572-2-RR



Basic Information

Catalog Number:

86572-2-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26028

GenBank Accession Number:

BC053544

GeneID (NCBI):

29801

UNIPROT ID:

Q9ULC8

Full Name:

zinc finger, DHHC-type containing 8

Calculated MW:

83 kDa

Observed MW:

70-85 kDa

Purification Method:

Protein A purification

CloneNo.:

251438A7

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : U2OS cells, HeLa cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86572-2-RR (ZDHHC8 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86572-2-RR against Human ZDHHC8 were performed. The affinity constant is 0.915 nM.