

For Research Use Only

# STX8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86580-2-PBS



## Basic Information

Catalog Number:

86580-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2850

GenBank Accession Number:

BC009713

GeneID (NCBI):

9482

UNIPROT ID:

Q9UNK0

Full Name:

syntaxin 8

Calculated MW:

236 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

251516H4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Syntaxin 8 concentration was determined to be 28.1 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 43.4 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02538-1, Syntaxin 8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86580-2-PBS. Detection antibody: 86580-1-PBS. Standard: Ag2850. Range: 0.469-60 ng/mL.