

For Research Use Only

STX8 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86580-3-PBS



Basic Information

Catalog Number:

86580-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2850

GenBank Accession Number:

BC009713

GeneID (NCBI):

9482

UNIPROT ID:

Q9UNK0

Full Name:

syntaxin 8

Calculated MW:

236 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

251516B7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Syntaxin 8 concentration was determined to be 46.48 ng/mL in A375 cell extract based on a 1.60 mg/mL extract load and 38.08 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02538-2, Human Syntaxin 8 Recombinant Matched Antibody Pair - PBS only. 86580-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2850. 86580-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.