

For Research Use Only

ORC6L Recombinant monoclonal antibody, PBS Only

Catalog Number:86591-1-PBS



Basic Information

Catalog Number:

86591-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11769

GenBank Accession Number:

BC063565

GeneID (NCBI):

23594

UNIPROT ID:

Q9Y5N6

Full Name:

origin recognition complex, subunit 6 like (yeast)

Calculated MW:

252 aa, 28 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

251479G7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86591-1-RR (ORC6L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86591-1-PBS in a different storage buffer formulation.



mouse testis tissue were subjected to SDS PAGE followed by western blot with 86591-1-RR (ORC6L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86591-1-PBS in a different storage buffer formulation.