

For Research Use Only

ORC6L Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86591-2-PBS



Basic Information

Catalog Number: 86591-2-PBS	GenBank Accession Number: BC063565	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 23594	CloneNo.: 251479G5
Isotype: IgG	UNIPROT ID: Q9Y5N6	
Immunogen Catalog Number: AG11769	Full Name: origin recognition complex, subunit 6 like (yeast)	
	Calculated MW: 252 aa, 28 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ORC6 concentration was determined to be 6.7 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load, 23.5 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02534-1, ORC6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86591-3-PBS. Detection antibody: 86591-2-PBS. Standard: Ag11769. Range: 0.312-40 ng/mL.