

For Research Use Only

TAX1BP1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 86602-3-PBS



Basic Information

Catalog Number: 86602-3-PBS	GenBank Accession Number: BC050358	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 8887	CloneNo.: 251529G8
Isotype: IgG	UNIPROT ID: Q86VP1	
Immunogen Catalog Number: AG5789	Full Name: Tax1 (human T-cell leukemia virus type I) binding protein 1	
	Calculated MW: 91 kDa	
	Observed MW: 90 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

TAX1BP1 (also named as T6BP and TXBP151) inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. It is degraded by caspase-3-like family proteins upon TNF-induced apoptosis. TAX1BP1 may also play a role in the pro-inflammatory cytokine IL-1 signaling cascade. TAX1BP1 binds TRAF6 but not TRAFs 1-5.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86602-3-RR (TAX1BP1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86602-3-PBS in a different storage buffer formulation.