

For Research Use Only

TAX1BP1 Recombinant monoclonal antibody

Catalog Number: 86602-3-RR



Basic Information

Catalog Number: 86602-3-RR	GenBank Accession Number: BC050358	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 8887	CloneNo.: 251529G8
Isotype: IgG	UNIPROT ID: Q86VP1	Recommended Dilutions: WB: 1:1000-1:6000
Immunogen Catalog Number: AG5789	Full Name: Tax1 (human T-cell leukemia virus type I) binding protein 1	
	Calculated MW: 91 kDa	
	Observed MW: 90 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells, A549 cells, U-937 cells
Species Specificity: human	

Background Information

TAX1BP1 (also named as T6BP and TXBP151) inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. It is degraded by caspase-3-like family proteins upon TNF-induced apoptosis. TAX1BP1 may also play a role in the pro-inflammatory cytokine IL-1 signaling cascade. TAX1BP1 binds TRAF6 but not TRAFs 1-5.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86602-3-RR (TAX1BP1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.