

For Research Use Only

ATXN2L Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86613-2-PBS



Basic Information

Catalog Number:

86613-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG20581

GenBank Accession Number:

BC082760

GeneID (NCBI):

11273

UNIPROT ID:

Q8WWM7

Full Name:

ataxin 2-like

Calculated MW:

1062 aa, 112 kDa

Purification Method:

Protein A purification

CloneNo.:

251578A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ATXN2L concentration was determined to be 102.6 ng/mL in A431 cell extract based on a 3.7 mg/mL extract load, 30.8 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load, 52.1 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02547-1, ATXN2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86613-1-PBS. Detection antibody: 86613-2-PBS. Standard: Ag20581. Range: 0.312-40 ng/mL.