

For Research Use Only

LRPPRC Recombinant antibody, PBS Only (Capture)

Catalog Number: 86616-1-PBS



Basic Information

Catalog Number: 86616-1-PBS	GenBank Accession Number: BC050311	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 10128	CloneNo.: 251577C1
Source: Rabbit	UNIPROT ID: P42704	
Isotype: IgG	Full Name: leucine-rich PPR-motif containing	
Immunogen Catalog Number: AG15452	Calculated MW: 1394 aa, 158 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LRPPRC concentration was determined to be 67.2 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load, 168.4 ng/mL in HCT 116 cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02559-1, LRPPRC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86616-1-PBS. Detection antibody: 86616-2-PBS. Standard: Ag15452. Range: 1.562-100 ng/mL.