

For Research Use Only

CD13 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86620-2-PBS



Basic Information

Catalog Number:

86620-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3748

GenBank Accession Number:

BC058928

GeneID (NCBI):

290

ENSEMBL Gene ID:

ENSG00000166825

UNIPROT ID:

P15144

Full Name:

alanyl (membrane) aminopeptidase

Calculated MW:

110 kDa

Purification Method:

Protein A purification

CloneNo.:

251501H7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD13 concentration was determined to be 726.6 ng/mL in human plasma.



Cytometric bead array standard curve of MP02566-1, CD13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86620-1-PBS. Detection antibody: 86620-2-PBS. Standard: Eg3748. Range: 0.625-80 ng/mL.