

For Research Use Only

INTS12 Recombinant antibody

Catalog Number: 86624-1-RR



Basic Information

Catalog Number:

86624-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9403

GenBank Accession Number:

BC014442

GeneID (NCBI):

57117

UNIPROT ID:

Q96CB8

Full Name:

integrator complex subunit 12

Calculated MW:

462 aa, 49 kDa

Observed MW:

49 kDa

Purification Method:

Protein A purification

CloneNo.:

251507A10

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, L02 cells

Background Information

INTS12 is a subunit of the Integrator complex that is thought to act as a reader to link the complex to chromatin.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86624-1-RR (INTS12 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.