

For Research Use Only

B4GALT1 Recombinant antibody, PBS Only (Capture)

Catalog Number:86640-1-PBS



Basic Information

Catalog Number:

86640-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001497

GeneID (NCBI):

2683

UNIPROT ID:

P15291

Full Name:

UDP-Gal:betaGlcNAc beta 1,4-galactosyltransferase, polypeptide 1

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

251527B5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00806-2, Mouse G-CSF Recombinant Matched Antibody Pair - PBS only. 86640-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0256. 86639-1-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL



Serum of sixteen mice was measured. The G-CSF concentration of detected samples was determined to be 594.3 pg/mL with a range of 189.4-987.2 pg/mL