

For Research Use Only

ROCK2(middle) Recombinant monoclonal antibody, PBS Only

Catalog Number:86649-4-PBS



Basic Information

Catalog Number:

86649-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16330

GenBank Accession Number:

BC111801

GeneID (NCBI):

9475

UNIPROT ID:

O75116

Full Name:

Rho-associated, coiled-coil containing protein kinase 2

Calculated MW:

77 kDa, 161 kDa

Observed MW:

161 kDa

Purification Method:

Protein A purification

CloneNo.:

251651F9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

ROCK2 is a serine/threonine kinase that regulates cytokinesis, smooth muscle contraction, the formation of actin stress fibers and focal adhesions, and the activation of the c-fos serum response element. It is an isozyme of ROCK1, is a target for the small GTPase Rho. ROCK2 is localized in both cytoplasm and nucleus, and the nucleus-localized ROCK2 can target p300 for phosphorylation to regulate its acetyltransferase activity (PMID:16574662). This antibody is specific to ROCK2.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86649-4-RR (ROCK2(middle) antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86649-4-PBS in a different storage buffer formulation.