

For Research Use Only

# MYH11 Recombinant antibody, PBS Only (Detector)

Catalog Number:86654-1-PBS



## Basic Information

<b>Catalog Number:</b> 86654-1-PBS	<b>GenBank Accession Number:</b> BC104906	<b>Purification Method:</b> Protein A purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 4629	<b>CloneNo.:</b> 251604D6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P35749	
<b>Isotype:</b> IgG	<b>Full Name:</b> myosin, heavy chain 11, smooth muscle	
<b>Immunogen Catalog Number:</b> AG16114	<b>Calculated MW:</b> 1979 aa, 228 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean SMMHC concentration was determined to be 15.6 ng/mL in mouse uterus tissue extract based on a 3.4 mg/mL extract load, 5.7 ng/mL in mouse eye tissue extract based on a 12.4 mg/mL extract load.



Cytometric bead array standard curve of MP02573-1, SMMHC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86654-2-PBS. Detection antibody: 86654-1-PBS. Standard: Ag16114. Range: 0.312-40 ng/mL.