

For Research Use Only

RHOB Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86663-1-PBS



Basic Information

Catalog Number:

86663-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5593

GenBank Accession Number:

BC066954

GeneID (NCBI):

388

UNIPROT ID:

P62745

Full Name:

ras homolog gene family, member B

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

251662F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RHOB concentration was determined to be 4.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02581-1, RHOB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86663-2-PBS. Detection antibody: 86663-1-PBS. Standard: Ag5593. Range: 0.781-100 ng/mL