

For Research Use Only

SAMD8 Recombinant monoclonal antibody

Catalog Number:86670-1-RR



Basic Information

Catalog Number:

86670-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25310

GenBank Accession Number:

BC146556

GeneID (NCBI):

142891

UNIPROT ID:

Q96LT4

Full Name:

sterile alpha motif domain containing 8

Observed MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

251533B11

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse heart tissue, rat heart tissue, human heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86670-1-RR (SAMD8 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.