

For Research Use Only

NHP2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86677-2-PBS



Basic Information

Catalog Number: 86677-2-PBS	GenBank Accession Number: BC000009	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 55651	CloneNo.: 251700E7
Isotype: IgG	UNIPROT ID: Q9NX24	
Immunogen Catalog Number: AG7251	Full Name: NHP2 ribonucleoprotein homolog (yeast)	
	Calculated MW: 17 kDa	

Applications

Tested Applications:
IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NHP2 concentration was determined to be 211.4 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02595-1, NHP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86677-3-PBS. Detection antibody: 86677-2-PBS. Standard: Ag7251. Range: 0.469-60 ng/mL.