

For Research Use Only

MAOB Recombinant monoclonal antibody

Catalog Number:86695-1-RR



Basic Information

Catalog Number: 86695-1-RR	GenBank Accession Number: BC022494	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 4129	CloneNo.: 251600B8
Isotype: IgG	UNIPROT ID: P27338	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:500-1:2000
Immunogen Catalog Number: AG3253	Full Name: monoamine oxidase B	
	Calculated MW: 520 aa, 59 kDa	
	Observed MW: 59 kDa	

Applications

Tested Applications:
WB, IHC, ELISA

Species Specificity:
human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : MCF-7 cells, mouse liver tissue, rat liver tissue

IHC : mouse liver tissue,

Background Information

The MAOB gene encodes a 520-amino acid protein with a molecular mass of 58.8 kD and shows 70% amino acid identity to MAOA. MAOA and MAOB are present in the outer mitochondrial membrane in the central nervous system and peripheral tissues. A comparison of highly purified human placental MAOA and human liver MAOB revealed that the A form of the enzyme is larger by 2 kDa and only one potential N-glycosylation site exists in each protein, with each site in a different relevant position (PMID:3387449).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




MCF-7 cells were subjected to SDS PAGE followed by western blot with 86695-1-RR (MAOB antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 86695-1-RR (MAOB antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 86695-1-RR (MAOB antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 86695-1-RR (MAOB antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).