

For Research Use Only

CXCL7/PPBP Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86698-2-PBS



Basic Information

Catalog Number:

86698-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3690

GenBank Accession Number:

NM_002704.3

GeneID (NCBI):

5473

UNIPROT ID:

P02775

Full Name:

pro-platelet basic protein (chemokine (C-X-C motif) ligand 7)

Calculated MW:

14kDa

Purification Method:

Protein A purification

CloneNo.:

250712B11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Plasma of eleven individual healthy human donors was measured. The human CXCL7/PPBP concentration of detected samples was determined to be 3,322.29 ng/mL with a range of 678.34 - 5,123.18 ng/mL



Sandwich ELISA standard curve of MP02754-1, Human CXCL7/PPBP Recombinant Matched Antibody Pair - PBS only. 86698-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3690. 86698-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL