

For Research Use Only

SFRS2 Recombinant monoclonal antibody, PBS Only

Catalog Number:86700-2-PBS



Basic Information

Catalog Number:

86700-2-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003016

GeneID (NCBI):

6427

UNIPROT ID:

Q01130

Full Name:

splicing factor, arginine/serine-rich 2

Calculated MW:

25 kDa

Observed MW:

25-35 kDa

Purification Method:

Protein A purification

CloneNo.:

251608H2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86700-2-RR (SFRS2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86700-2-PBS in a different storage buffer formulation.